

Flowing Gas

Current State of Detector: Detector is being pumped

- Valves that are closed: $0a_{bot}$ (bottle)
- Valves that are open: $0a_{reg}$ (regulator), 1a, 3, 4, 5, 6, 7, 8, 9, 10

Flowing gas:

- Close Valves: $0a_{bot}$ (should already be closed), $0a_{reg}$, 1a, 5, 8, 9, 10
 - $0a_{bot}$, $0a_{reg}$ closes the gas bottle and regulator
 - 1a closes the gas line
 - 5 closes the detector return
 - Eventually, we will slowly/barely open 5 to get the correct flow rate
 - 6, 7 control flow through source holder
 - 8 closes one of the bypasses
 - 9 closes one of the bypasses
 - 10 closes the large pump-out
 - This should cause pressure to increase slightly
- Open Valves: 3, 4, 6, 7 (all should already be open)
 - 3 is the main line
 - 4 is the detector source
 - With 4 open and 5 closed we can let gas into the detector and control flow rate by finely adjusting 5
- Open $0a_{bot}$
- Open $0a_{reg}$
 - Set to ~5 psi
 - Make sure blue release valve is open
- Open valve 1a
- Set mass-flow controller to desired pressure
 - Typical operating pressure is 800 Torr
- Adjust valve 5 to desired flow rate
 - Flow rate should be ~20 sccm (standard cubic centimeters per minute) for normal operation
 - Flow rate should be ~1000 sccm if using the GHS source holder

Pumping

Current State of Detector: Gas is Flowing through the Detector

- Valves that are closed: 8, 9, 10
- Valves that are open: $0a_{bot}$, $0a_{reg}$, 1a, 3, 4, 5 (partially), 6, 7

Pumping:

- Set mass-flow controller to zero pressure
- Close Valves: $0a_{bot}$, 4, 5, (9, 10 should already be closed)
 - $0a_{bot}$ stops gas
 - 4, 5 seals off the detector
- Open Valves: 8, ($0a_{reg}$, 1a, 3, 6, 7 should already be open)
- Open valve 9
- Open valve 5 for slow pumping
 - Don't want to create too much wind in the detector, which is why 4 is closed and we are only slowly opening 5.
- When pressure in chamber is below 200 Torr, open valves 4, 10, and 5 all the way